



MZT2A mouse mAb

Catalog No	YP-mAb-09844
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MZT2A FAM128A MOZART2A
Immunogen	Synthesized peptide derived from human MZT2A AA range: 72-122
Specificity	This antibody detects endogenous levels of human MZT2A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AMS2;EC:2.5.1.6;MAT 2;MATA2;MATII;MZT2A;MZT2A FAM128A MOZART2A; Mitotic-spindle organizing protein 2A (Mitotic-spindle organizing protein associated with a ring of gamma-tubulin 2A)
Observed Band	47 kDa, 49 kDa
Calculated Molecular Weight	44 kDa
Cell Pathway	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome. Cytoplasm, cytoskeleton, spindle.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images