



C170L mouse mAb

Catalog No	YP-mAb-09847
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CEP170P1 CEP170L KIAA0470L
Immunogen	Synthesized peptide derived from human C170L AA range: 13-63
Specificity	This antibody detects endogenous levels of human C170L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C170L;CEP170P1 CEP170L KIAA0470L;Cep170-like protein (CEP170 pseudogene 1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images