



PHYIP mouse mAb

Catalog No	YP-mAb-09849
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PHYHIP DYRK1AP3 KIAA0273
Immunogen	Synthesized peptide derived from human PHYIP AA range: 1-51
Specificity	This antibody detects endogenous levels of human PHYIP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DYRK1AP3;KIAA0273;PAHX AP;PAHX AP1;PAHXAP1;PHYHIP;PHYIP; PHYHIP DYRK1AP3 KIAA0273;Phytanoyl-CoA hydroxylase-interacting protein (Phytanoyl-CoA hydroxylase-associated protein 1) (PAHX-AP1) (PAHXAP1)
Observed Band	35 kDa, 43-45 kDa
Calculated Molecular Weight	330 aa, 38 kDa
Cell Pathway	
Tissue Specificity	Highly expressed in the brain.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images