



KIFA3 mouse mAb

Catalog No	YP-mAb-09850
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KIFAP3 KIF3AP SMAP
Immunogen	Synthesized peptide derived from human KIFA3 AA range: 65-115
Specificity	This antibody detects endogenous levels of human KIFA3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Kinesin-associated protein 3;kinesin associated protein 3;KIF3AP;KAP-3;KAP3;KIFA3;KIFAP3 KIF3AP SMAP;Kinesin-associated protein 3 (KAP-3) (KAP3) (Smg GDS-associated protein)
Observed Band	91-100 kDa
Calculated Molecular Weight	792 aa, 91 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images