



MBD5 mouse mAb

Catalog No	YP-mAb-09851
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MBD5 KIAA1461
Immunogen	Synthesized peptide derived from human MBD5 AA range: 1057-1107
Specificity	This antibody detects endogenous levels of human MBD5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA1461;MBD5;MRD1;MBD5 KIAA1461;Methyl-CpG-binding domain protein 5 (Methyl-CpG-binding protein MBD5)
Observed Band	160 kDa 160 kDa
Calculated Molecular Weight	292 aa, 31 kDa,
Cell Pathway	[Isoform 1]: Nucleus. Chromosome. Associated with pericentric heterochromatin.; [Isoform 2]: Nucleus. Not associated with pericentric heterochromatin.
Tissue Specificity	Detected in heart, placenta, liver, skeletal muscle, kidney and pancreas.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images