



EDAD mouse mAb

Catalog No	YP-mAb-09852
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	EDARADD
Immunogen	Synthesized peptide derived from human EDAD AA range: 108-158
Specificity	This antibody detects endogenous levels of human EDAD
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EDAD;EDARADD;Ectodysplasin-A receptor-associated adapter protein (EDAR-associated death domain protein) (Protein crinkled homolog)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	Detected in adult pancreas, placenta and fetal skin, and at lower levels in lung, thymus, prostate and testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images