



U2AFL mouse mAb

Catalog No	YP-mAb-09854
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZRSR1 U2AF1-RS1 U2AF1L1 U2AF1P U2AF1RS1 U2AFBPL
Immunogen	Synthesized peptide derived from human U2AFL AA range: 131-181
Specificity	This antibody detects endogenous levels of human U2AFL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	U2AFL;ZRSR1 U2AF1-RS1 U2AF1L1 U2AF1P U2AF1RS1 U2AFBPL;U2 small nuclear ribonucleoprotein auxiliary factor 35 kDa subunit-related protein 1 (CCCH type zinc finger;RNA-binding motif and serine/arginine rich protein 1) (U2(RNU2) small nuclear RNA auxiliary; RNA-binding motif and serine/arginine rich protein 1) (U2(RNU2) small nuclear RNA auxiliary
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images