



PKHH2 mouse mAb

Catalog No	YP-mAb-09856
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PLEKHH2 KIAA2028
Immunogen	Synthesized peptide derived from human PKHH2 AA range: 565-615
Specificity	This antibody detects endogenous levels of human PKHH2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PKHH2;PLEKHH2 KIAA2028;Pleckstrin homology domain-containing family H member 2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Cytoplasm, cytoskeleton . Cell membrane ; Peripheral membrane protein ; Cytoplasmic side . Cell projection, lamellipodium . Localizes to foot process of podocytes. Localization to peripheral regions of lamellipodia seems to be dependent on PI3K.
Tissue Specificity	Kidney. Reduced expression in patients with focal segmental glomerulosclerosis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images