



PLB1 mouse mAb

Catalog No	YP-mAb-09857
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PLB1 PLB
Immunogen	Synthesized peptide derived from human PLB1 AA range: 72-122
Specificity	This antibody detects endogenous levels of human PLB1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PLB;PLB/LIP;PLB1;PLB1 PLB;Phospholipase B1;membrane-associated (Phospholipase B) (hPLB) (Phospholipase B/lipase) (PLB/LIP) [Includes: Phospholipase A2 (EC 3.1.1.4);Lysophospholipase (EC 3.1.1.5)]; membrane-associated (Phospholipase B) (hPLB) (Phospholipase B/lipase) (PLB/LIP) [Includes: Phospholipase A2 (EC 3.1.1.4)]
Observed Band	170kDa
Calculated Molecular Weight	163kDa
Cell Pathway	Apical cell membrane ; Single-pass type I membrane protein . Present in the intestinal brush border membranes. .
Tissue Specificity	Expressed in the epidermis (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images