



ISK1 mouse mAb

Catalog No	YP-mAb-09859
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SPINK1 PSTI
Immunogen	Synthesized peptide derived from human ISK1 AA range: 26-75
Specificity	This antibody detects endogenous levels of human ISK1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ISK1;SPINK1 PSTI;Pancreatic secretory trypsin inhibitor (Serine protease inhibitor Kazal-type 1) (Tumor-associated trypsin inhibitor) (TATI)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images