



D2HDH mouse mAb

Catalog No	YP-mAb-09865
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	D2HGDH D2HGD
Immunogen	Synthesized peptide derived from human D2HDH AA range: 199-249
Specificity	This antibody detects endogenous levels of human D2HDH
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	D2HDH;D2HGDH D2HGD;D-2-hydroxyglutarate dehydrogenase;mitochondrial (EC 1.1.99.-); mitochondrial (EC 1.1.99.-)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images