



PARVG mouse mAb

Catalog No	YP-mAb-09869
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PARVG
Immunogen	Synthesized peptide derived from human PARVG AA range: 113-163
Specificity	This antibody detects endogenous levels of human PARVG
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	parvin;gamma;Gamma-parvin;Gamma parvin;PARVG
Observed Band	37 kDa
Calculated Molecular Weight	331 aa, 37 kDa
Cell Pathway	Cell junction, focal adhesion. Cell membrane ; Peripheral membrane protein ; Cytoplasmic side . Cytoplasm, cytoskeleton . Constituent of focal adhesions. .
Tissue Specificity	Expressed predominantly in lymphoid organs, including spleen, thymus, lymph node, bone marrow and peripheral blood leukocytes and moderately in the digestive tract, including stomach, duodenum, jejunum, ileum, ileocecum and appendix, as well as in lung and liver. Also expressed in tumors, but at a lower level than in the corresponding normal tissues.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images