



KDM4E mouse mAb

Catalog No	YP-mAb-09871
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	KDM4E KDM4DL
Immunogen	Synthesized peptide derived from human KDM4E AA range: 338-388
Specificity	This antibody detects endogenous levels of human KDM4E
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KDM4E;KDM4E KDM4DL;Lysine-specific demethylase 4E (EC 1.14.11.-) (KDM4D-like protein) (Lysine-specific demethylase 4D-like)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images