



NUD19 mouse mAb

Catalog No	YP-mAb-09872
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NUDT19
Immunogen	Synthesized peptide derived from human NUD19 AA range: 196-246
Specificity	This antibody detects endogenous levels of human NUD19
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NUD19;NUDT19;Nucleoside diphosphate-linked moiety X motif 19;mitochondrial (Nudix motif 19) (EC 3.6.1.-); mitochondrial (Nudix motif 19) (EC 3.6.1.-)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Peroxisome .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images