



GAK3 mouse mAb

Catalog No	YP-mAb-09876
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	Synthesized peptide derived from human GAK3 AA range: 137-187
Specificity	This antibody detects endogenous levels of human GAK3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GAK3;HERV-K_19q12 provirus ancestral Gag polyprotein (Gag polyprotein) (HERV-K(C19) Gag protein) [Cleaved into: Matrix protein;Capsid protein; Nucleocapsid protein]
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane; Lipid-anchor. Cytoplasmic membrane (in a transfection system).
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images