



BCAS1 mouse mAb

Catalog No	YP-mAb-09877
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	BCAS1 AIBC1 NABC1
Immunogen	Synthesized peptide derived from human BCAS1 AA range: 262-312
Specificity	This antibody detects endogenous levels of human BCAS1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Breast carcinoma-amplified sequence 1;NABC1;BCAS1;BCAS1 AIBC1 NABC1; Breast carcinoma-amplified sequence 1 (Amplified and overexpressed in breast cancer) (Novel amplified in breast cancer 1)
Observed Band	68 kDa
Calculated Molecular Weight	584 aa, 62 kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	Highly expressed in the brain and, more specifically, in oligodendrocytes (at protein level). Expressed in the prostate, and at lower levels in testis, intestine and colon. Overexpressed in most breast cancer cell lines and down-regulated in some colorectal tumors.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images