



LSM3 mouse mAb

Catalog No	YP-mAb-09878
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LSM3 MDS017
Immunogen	Synthesized peptide derived from human LSM3 AA range: 25-75
Specificity	This antibody detects endogenous levels of human LSM3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	YLR438C;USS2;U6 snRNA-associated Sm-like protein LSm3;SMX4;MDS017; LSM3;LSM3 MDS017
Observed Band	12 kDa
Calculated Molecular Weight	102 aa, 12 kDa
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images