



MKLN1 mouse mAb

Catalog No	YP-mAb-09882
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MKLN1
Immunogen	Synthesized peptide derived from human MKLN1 AA range: 61-111
Specificity	This antibody detects endogenous levels of human MKLN1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TWA2;MKLN1;Muskelin
Observed Band	100kDa
Calculated Molecular Weight	85kDa
Cell Pathway	Cytoplasm . Cytoplasm, cytosol . Nucleus, nucleoplasm . Cell projection, ruffle . Cytoplasm, cell cortex . Cell junction, synapse . Cell junction, synapse, postsynapse . Colocalizes with GABRA1 at synapses and in postsynaptic regions. Colocalizes with actin fibers in the cell cortex. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images