



LMAN2 mouse mAb

Catalog No	YP-mAb-09884
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LMAN2 C5orf8
Immunogen	Synthesized peptide derived from human LMAN2 AA range: 227-277
Specificity	This antibody detects endogenous levels of human LMAN2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C5orf8;Glycoprotein GP36b;Lectin mannose-binding 2;Vesicular integral-membrane protein 36;Vesicular integral-membrane protein VIP36;LMAN2; LMAN2 C5orf8;Vesicular integral-membrane protein VIP36 (Glycoprotein GP36b) (Lectin mannose-binding 2) (Vesicular integral-membrane protein 36) (VIP36)
Observed Band	35-36 kDa
Calculated Molecular Weight	356 aa, 40 kDa
Cell Pathway	Endoplasmic reticulum-Golgi intermediate compartment membrane ; Single-pass type I membrane protein . Golgi apparatus membrane ; Single-pass membrane protein . Endoplasmic reticulum membrane ; Single-pass type I membrane protein .
Tissue Specificity	Ubiquitous.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images