



RRMJ3 mouse mAb

Catalog No	YP-mAb-09885
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FTSJ3 SB92
Immunogen	Synthesized peptide derived from human RRMJ3 AA range: 628-678
Specificity	This antibody detects endogenous levels of human RRMJ3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RRMJ3;FTSJ3 SB92;pre-rRNA processing protein FTSJ3 (EC 2.1.1.-) (Protein ftsJ homolog 3) (Putative rRNA methyltransferase 3) (rRNA (uridine-2'-O)-methyltransferase 3)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images