



CD022 mouse mAb

Catalog No	YP-mAb-09886
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C4orf22
Immunogen	Synthesized peptide derived from human CD022 AA range: 153-203
Specificity	This antibody detects endogenous levels of human CD022
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CD022;C4orf22;Uncharacterized protein C4orf22
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Nucleus . Mainly cytoplasmic. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images