



SO1B7 mouse mAb

Catalog No	YP-mAb-09887
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SLCO1B7 LST3 LST3TM12
Immunogen	Synthesized peptide derived from human SO1B7 AA range: 270-320
Specificity	This antibody detects endogenous levels of human SO1B7
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SO1B7;SLCO1B7 LST3 LST3TM12;Putative solute carrier organic anion transporter family member 1B7 (Liver-specific organic anion transporter 3) (LST-3)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images