



ARK72 mouse mAb

Catalog No	YP-mAb-09890
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	AKR7A2 AFAR AFAR1 AKR7
Immunogen	Synthesized peptide derived from human ARK72 AA range: 219-269
Specificity	This antibody detects endogenous levels of human ARK72
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ARK72;AKR7A2 AFAR AFAR1 AKR7;Aflatoxin B1 aldehyde reductase member 2 (EC 1.1.1.n11) (AFB1 aldehyde reductase 1) (AFB1-AR 1) (Aldoketoreductase 7) (Succinic semialdehyde reductase) (SSA reductase)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus . Cytoplasm .
Tissue Specificity	Detected in brain, liver, small intestine and testis, and at lower levels in heart, prostate, skeletal muscle and spleen. Detected in kidney proximal and distal tubules, endothelial cells lining the Bowman's capsules and some cysts. Detected at low levels in lung and pancreas (at protein level). Widely expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images