



TXND9 mouse mAb

Catalog No	YP-mAb-09893
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TXNDC9 APACD
Immunogen	Synthesized peptide derived from human TXND9 AA range: 149-199
Specificity	This antibody detects endogenous levels of human TXND9
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	APACD;ATP-binding protein associated with cell differentiation;PHLP3;Protein 1 4;Protein 1-4;TXND9;TXNDC9 APACD;Thioredoxin domain-containing protein 9 (ATP-binding protein associated with cell differentiation) (Protein 1-4)
Observed Band	25-27 kDa
Calculated Molecular Weight	226 aa, 27 kDa
Cell Pathway	Cytoplasm . Nucleus . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Midbody . Co-localizes with beta-tubulin in the centrosome. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images