



DOCK9 mouse mAb

Catalog No	YP-mAb-09894
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DOCK9 KIAA1058
Immunogen	Synthesized peptide derived from human DOCK9 AA range: 1516-1566
Specificity	This antibody detects endogenous levels of human DOCK9
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Zizimin-1;ZIZ1;dedicator of cytokinesis 9;Cdc42 guanine nucleotide exchange factor zizimin-1;DOCK9;DOCK9 KIAA1058;Dedicator of cytokinesis protein 9 (Cdc42 guanine nucleotide exchange factor zizimin-1)
Observed Band	200-236 kDa
Calculated Molecular Weight	239 kDa
Cell Pathway	Endomembrane system . Associated with membranes. .
Tissue Specificity	Widely expressed, with highest expression in heart and placenta. Expressed at intermediate level in kidney, brain, lung and skeletal muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images