



MUC19 mouse mAb

Catalog No	YP-mAb-09900
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MUC19
Immunogen	Synthesized peptide derived from human MUC19 AA range: 64-114
Specificity	This antibody detects endogenous levels of human MUC19
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MUC19;Mucin-19 (MUC-19)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	Expressed corneal epithelial cells, conjunctival goblet and epithelial cells and lacrimal gland cells (at protein level). Expressed by mucous cells of the submandibular gland and submucosal gland of the trachea. Expressed by middle ear epithelial cells.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images