



PLGB mouse mAb

Catalog No	YP-mAb-09901
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PLGLB1 PLGL PRGB; PLGLB2 PLGP1
Immunogen	Synthesized peptide derived from human PLGB AA range: 14-64
Specificity	This antibody detects endogenous levels of human PLGB
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GPI mannosyltransferase 3;GPI mannosyltransferase III;GPI MT III;PIG B;PIGB; PLGB;PLGLB1 PLGL PRGB;PLGLB2 PLGP1;Plasminogen-like protein B (Plasminogen-related protein B)
Observed Band	50-55 kDa
Calculated Molecular Weight	65 kDa
Cell Pathway	Secreted.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images