



ZN263 mouse mAb

Catalog No	YP-mAb-09903
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF263 FPM315 ZKSCAN12
Immunogen	Synthesized peptide derived from human ZN263 AA range: 557-607
Specificity	This antibody detects endogenous levels of human ZN263
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZKSCAN12;Zinc finger protein FPM315;Zinc finger protein 263;ZFP263;FPM315; ZN263;ZNF263 FPM315 ZKSCAN12;Zinc finger protein 263 (Zinc finger protein FPM315) (Zinc finger protein with KRAB and SCAN domains 12)
Observed Band	77 kDa
Calculated Molecular Weight	683 aa, 77 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in heart, brain, placenta, lung, liver, skeletal muscle, kidney, pancreas, spleen, thymus, prostate, testis, ovary, small intestine, colon and leukocyte.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images