



MACD2 mouse mAb

Catalog No	YP-mAb-09908
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MACROD2 C20orf133
Immunogen	Synthesized peptide derived from human MACD2 AA range: 166-216
Specificity	This antibody detects endogenous levels of human MACD2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MACD2;MACROD2 C20orf133;O-acetyl-ADP-ribose deacetylase MACROD2 (EC 3.5.1.-) (MACRO domain-containing protein 2)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Recruited to DNA lesions, probably via mono-APD-ribosylated proteins. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images