



QPCTL mouse mAb

Catalog No	YP-mAb-09911
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	QPCTL
Immunogen	Synthesized peptide derived from human QPCTL AA range: 152-202
Specificity	This antibody detects endogenous levels of human QPCTL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Glutaminyl cyclase like;Glutaminyl peptide cyclotransferase like;Glutaminyl peptide cyclotransferase like protein;Glutaminyl-peptide cyclotransferase-like protein;Golgi-resident glutaminyl-peptide cyclotransferase;gQC;isoQC;QPCTL; QPCTL_HUMAN;Glutaminyl-peptide cyclotransferase-like protein (EC 2.3.2.5) (Golgi-resident glutaminyl-peptide cyclotransferase) (isoQC) (gQC)
Observed Band	42kD
Calculated Molecular Weight	43kD(Calculated).
Cell Pathway	Golgi apparatus membrane ; Single-pass type I membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images