



NFIA mouse mAb

Catalog No	YP-mAb-09919
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NFIA KIAA1439
Immunogen	Synthesized peptide derived from human NFIA AA range: 391-441
Specificity	This antibody detects endogenous levels of human NFIA
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CTF;KIAA1439;NF I/A;NF1 A;NFI A;NFI L;NFIA;Nuclear factor 1 A type;Nuclear factor 1/A;nuclear factor I/A;TGGCA binding protein;NFIA KIAA1439;Nuclear factor 1 A-type (NF1-A) (Nuclear factor 1/A) (CCAAT-box-binding transcription factor) (CTF) (Nuclear factor I/A) (NF-I/A) (NFI-A) (TGGCA-binding protein)
Observed Band	60-70 kDa
Calculated Molecular Weight	498 aa, 55 kDa
Cell Pathway	Nucleus.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images