



ASB14 mouse mAb

Catalog No	YP-mAb-09921
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ASB14
Immunogen	Synthesized peptide derived from human ASB14 AA range: 507-557
Specificity	This antibody detects endogenous levels of human ASB14
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ASB 14;ASB14;Ankyrin repeat and SOCS box protein 14 (ASB-14)
Observed Band	65 kDa
Calculated Molecular Weight	587 aa, 65 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images