



SHSA9 mouse mAb

Catalog No	YP-mAb-09922
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SHISA9
Immunogen	Synthesized peptide derived from human SHSA9 AA range: 352-402
Specificity	This antibody detects endogenous levels of human SHSA9
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SHSA9;SHISA9;Protein shisa-9
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell projection, dendritic spine membrane ; Single-pass type I membrane protein . Cell junction, synapse .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images