



HXD4 mouse mAb

Catalog No	YP-mAb-09923
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HOXD4 HOX4B
Immunogen	Synthesized peptide derived from human HXD4 AA range: 99-149
Specificity	This antibody detects endogenous levels of human HXD4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HXD4;HOXD4 HOX4B;Homeobox protein Hox-D4 (Homeobox protein HHO.C13) (Homeobox protein Hox-4B) (Homeobox protein Hox-5.1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images