



WDR41 mouse mAb

Catalog No	YP-mAb-09925
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	WDR41 MSTP048
Immunogen	Synthesized peptide derived from human WDR41 AA range: 200-250
Specificity	This antibody detects endogenous levels of human WDR41
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	WD repeat-containing protein 41;WD repeat domain 41;MSTP048;WDR41;WDR41 MSTP048
Observed Band	52 kDa
Calculated Molecular Weight	459 aa, 52 kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images