



RNZ1 mouse mAb

Catalog No	YP-mAb-09926
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ELAC1 D29
Immunogen	Synthesized peptide derived from human RNZ1 AA range: 126-176
Specificity	This antibody detects endogenous levels of human RNZ1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RNZ1;ELAC1 D29;Zinc phosphodiesterase ELAC protein 1 (EC 3.1.26.11) (Deleted in Ma29) (ElaC homolog protein 1) (Ribonuclease Z 1) (RNase Z 1) (tRNA 3 endonuclease 1) (tRNase Z 1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytosol . Nucleus . Mainly cytosolic.
Tissue Specificity	Widely expressed. Expressed in heart, brain, placenta, lung, liver, skeletal muscle, kidney and pancreas.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images