



ADAP2 mouse mAb

Catalog No	YP-mAb-09928
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ADAP2 CENTA2
Immunogen	Synthesized peptide derived from human ADAP2 AA range: 301-351
Specificity	This antibody detects endogenous levels of human ADAP2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ADAP2;ArfGAP with dual PH domains 2;cent b;CENTA2;Centaurin alpha 2;Cnt a2;HSA272195;ADAP2 CENTA2;Arf-GAP with dual PH domain-containing protein 2 (Centaurin-alpha-2) (Cnt-a2)
Observed Band	44 kDa
Calculated Molecular Weight	380 aa, 44 kDa
Cell Pathway	Cytoplasm . Cell membrane . Constitutively associated with the plasma membrane. Excluded from the nucleus.
Tissue Specificity	Highly expressed in placenta, spleen, kidney, skeletal muscle and adrenal gland. Weakly expressed in thyroid, liver, heart, lung, small intestine, peripheral blood leukocytes. Not detected in spinal cord, brain, stomach, trachea, colon, lymph node and bone marrow.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images