



SKAP mouse mAb

Catalog No	YP-mAb-09930
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SKAP C15orf23 TRAF4AF1 HSD11
Immunogen	Synthesized peptide derived from human SKAP AA range: 34-84
Specificity	This antibody detects endogenous levels of human SKAP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SKAP;SKAP C15orf23 TRAF4AF1 HSD11;Small kinetochore-associated protein (SKAP) (Kinetochore-localized astrin-binding protein) (Kinastrin) (TRAF4-associated factor 1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Chromosome, centromere, kinetochore . Cytoplasm, cytoskeleton, spindle pole . Cytoplasm, cytoskeleton, microtubule organizing center . Colocalizes with microtubules around centrosomes in prophase and with the mitotic spindle at prometaphase and metaphase. From late prometaphase to anaphase, is highly concentrated on kinetochores. Located at the kinetochore-microtubule interface. The astrin (SPAG5)-kinastrin (SKAP) complex localizes to the microtubule plus ends (PubMed:23035123). .
Tissue Specificity	Widely expressed, including in skin.
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images