



CNIH3 mouse mAb

Catalog No	YP-mAb-09931
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CNIH3
Immunogen	Synthesized peptide derived from human CNIH3 AA range: 44-94
Specificity	This antibody detects endogenous levels of human CNIH3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CNIH3;Protein cornichon homolog 3
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell junction, synapse, postsynaptic cell membrane ; Multi-pass membrane protein . Also localizes to the cell membrane of extrasynaptic sites (dendritic shafts, spines of pyramidal cells). .
Tissue Specificity	Expression is up-regulated in dorsolateral prefrontal cortex of patients with schizophrenia (postmortem brain study).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images