



ARH mouse mAb

Catalog No	YP-mAb-09935
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LDLRAP1 ARH
Immunogen	Synthesized peptide derived from human ARH AA range: 236-286
Specificity	This antibody detects endogenous levels of human ARH
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LDLRAP1;ARH1;ARH2;FHCB1;FHCB2;ARH;LDLRAP1 ARH;Low density lipoprotein receptor adapter protein 1 (Autosomal recessive hypercholesterolemia protein)
Observed Band	34 kDa
Calculated Molecular Weight	308 aa, 34 kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	Expressed at high levels in the kidney, liver, and placenta, with lower levels detectable in brain, heart, muscle, colon, spleen, intestine, lung, and leukocytes.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images