



PERQ1 mouse mAb

Catalog No	YP-mAb-09937
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GIGYF1 CDS2 PERQ1 PP3360
Immunogen	Synthesized peptide derived from human PERQ1 AA range: 684-734
Specificity	This antibody detects endogenous levels of human PERQ1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PERQ1;GIGYF1 CDS2 PERQ1 PP3360;PERQ amino acid-rich with GYF domain-containing protein 1 (GRB10-interacting GYF protein 1)
Observed Band	114kD
Calculated Molecular Weight	114kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images