



ALKB6 mouse mAb

Catalog No	YP-mAb-09938
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ALKBH6 ABH6
Immunogen	Synthesized peptide derived from human ALKB6 AA range: 17-67
Specificity	This antibody detects endogenous levels of human ALKB6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ALKB6;ALKBH6 ABH6;Probable alpha-ketoglutarate-dependent dioxygenase ABH6 (EC 1.14.11.-) (Alkylated DNA repair protein alkB homolog 6)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	Widely expressed, with highest expression in testis and pancreas.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images