



LRC58 mouse mAb

Catalog No	YP-mAb-09939
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LRRC58
Immunogen	Synthesized peptide derived from human LRC58 AA range: 292-342
Specificity	This antibody detects endogenous levels of human LRC58
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LRC58;LRRC58;Leucine-rich repeat-containing protein 58
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images