



SPEE mouse mAb

Catalog No	YP-mAb-09942
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SRM SPS1 SRML1
Immunogen	Synthesized peptide derived from human SPEE AA range: 235-285
Specificity	This antibody detects endogenous levels of human SPEE
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SPEE;SRM SPS1 SRML1;Spermidine synthase (SPDSY) (EC 2.5.1.16) (Putrescine aminopropyltransferase)
Observed Band	33kD
Calculated Molecular Weight	33kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images