



CV015 mouse mAb

Catalog No	YP-mAb-09945
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C22orf15
Immunogen	Synthesized peptide derived from human CV015 AA range: 59-109
Specificity	This antibody detects endogenous levels of human CV015
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CV015;C22orf15;Uncharacterized protein C22orf15 (Protein N27C7-3)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images