



LCE3B mouse mAb

Catalog No	YP-mAb-09948
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LCE3B LEP14
Immunogen	Synthesized peptide derived from human LCE3B AA range: 8-58
Specificity	This antibody detects endogenous levels of human LCE3B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LCE3B;LCE3B LEP14;Late cornified envelope protein 3B (Late envelope protein 14)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Skin-specific. Expression was readily detected in adult trunk skin, adult arm skin, fetal skin, penal skin, vulva, esophagus and tongue. Not expressed in the cervix, rectum, lung, colon, or placenta.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images