



CCD89 mouse mAb

Catalog No	YP-mAb-09951
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CCDC89 BOIP
Immunogen	Synthesized peptide derived from human CCD89 AA range: 293-343
Specificity	This antibody detects endogenous levels of human CCD89
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CCD89;CCDC89 BOIP;Coiled-coil domain-containing protein 89 (Bc8 orange-interacting protein)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Nucleus . Uniformly distributed within the cell, but becomes recruited to the nucleus upon binding to HEY1. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images