



TMHS mouse mAb

Catalog No	YP-mAb-09952
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	LHFPL5 TMHS
Immunogen	Synthesized peptide derived from human TMHS AA range: 28-78
Specificity	This antibody detects endogenous levels of human TMHS
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TMHS;LHFPL5 TMHS;Tetraspan membrane protein of hair cell stereocilia (Lipoma HMGIC fusion partner-like 5 protein)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Multi-pass membrane protein . Efficient localization to the plasma membrane requires the presence of PCDH15. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images