



LMBL3 mouse mAb

Catalog No	YP-mAb-09955
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	L3MBTL3 KIAA1798 MBT1
Immunogen	Synthesized peptide derived from human LMBL3 AA range: 128-178
Specificity	This antibody detects endogenous levels of human LMBL3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LMBL3;L3MBTL3 KIAA1798 MBT1;Lethal(3)malignant brain tumor-like protein 3 (H-I(3)mbt-like protein 3) (L(3)mbt-like protein 3) (MBT-1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images