



WF10B mouse mAb

Catalog No	YP-mAb-09956
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	WFDC10B WAP12
Immunogen	Synthesized peptide derived from human WF10B AA range: 14-64
Specificity	This antibody detects endogenous levels of human WF10B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	WF10B;WFDC10B WAP12;Protein WFDC10B
Observed Band	
Calculated Molecular Weight	8kD
Cell Pathway	Secreted .
Tissue Specificity	Ubiquitously expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images